

Safety Instruction for Raman unit.

Read carefully before start operating with device!!!

Laser is ON if power supply turned on and Safe key is in the ON position.

Raman unit contains laser with output radiation at 532nm and 250mW. It can be dangerous for eyes!!!

Raman unit comes with safety goggles.

- **Never open Raman unit cover if laser is ON.**
- **Safety goggles should be put on when the laser is ON**
- **Everybody around (in direct visibility to the Raman unit) should wear safety goggles if laser is ON.**
- **Never work with Raman unit if cover of Raman unit or other box panels are broken or safety goggles broken.**

You can purchase extra safety goggles if you need several persons operate with Raman unit simultaneously